

Anti-Phospho-WASP (Y290) Rabbit Polyclonal Antibody



Catalog #: U0597

Aliases

IMD2; Wiskott-Aldrich syndrome protein; WASp

Background

Gene Name: WAS

NCBI Gene Entry: [7454](#)

UniProt Entry: [P42768](#)

Application Information

Molecular Weight: Predicted, 52 kDa; observed, 66,60 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB), immunohistochemistry (IHC), immunocytochemistry (IC), flow cytometry (FCM)

Immunogen

A synthesized peptide derived from human Phospho-WASP (Y290)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Immunohistochemistry (IHC): 1:100-1:200

Immunocytochemistry (IC): 1:100-1:500

Flow Cytometry (FCM) 1:100-1:200

Note: This product is for research use only.

Validation Data




SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX:+1-540-855-7041

WWW.GENUINBIOTECH.COM