

Anti-Phospho-Aconitase 1 (S711) Rabbit Polyclonal Antibody



Catalog #: U0600

Aliases

IREB1; Cytoplasmic aconitate hydratase; Aconitase; Citrate hydro-lyase; Ferritin repressor protein; Iron regulatory protein 1; IRP1; Iron-responsive element-binding protein 1; IRE-BP 1

Background

Gene Name: ACO1
NCBI Gene Entry: 48
UniProt Entry: P21399

Application Information

Molecular Weight: Predicted, 98 kDa; observed, 100 kDa
Clonality: Rabbit polyclonal antibody
Species Reactivity: Human, mouse, rat, bovine, chicken, rabbit, zebrafish
Applications Tested: Western blotting (WB), immunohistochemistry (IHC)

Immunogen

A synthesized peptide derived from human Phospho-Aconitase 1 (S711)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000
Immunohistochemistry (IHC): 1:100-1:200

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041



ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-Phospho-Aconitase 1 (S711) Rabbit Polyclonal Antibody

PAGE 2



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM