

Anti-Phospho-RhoGDI alpha (S174) Rabbit Polyclonal Antibody



Catalog #: U0606

Aliases

GDIA1; Rho GDP-dissociation inhibitor 1; Rho GDI 1; Rho-GDI alpha

Background

Gene Name: ARHGDIA

NCBI Gene Entry: [396](#)

UniProt Entry: [P52565](#)

Application Information

Molecular Weight: Predicted, 23 kDa; observed, 26 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, monkey

Applications Tested: Western blotting (WB), immunohistochemistry (IHC)

Immunogen

A synthesized peptide derived from human Phospho-RhoGDI alpha (S174)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Immunohistochemistry (IHC): 1:100-1:200

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041



ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-Phospho-RhoGDI alpha (S174) Rabbit Polyclonal Antibody

PAGE 2



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM