

Anti-Phospho-BCR (Y360) Rabbit Polyclonal Antibody



Catalog #: U0610

Aliases

BCR1; D22S11; Breakpoint cluster region protein; Renal carcinoma antigen NY-REN-26

Background

Gene Name: BCR

NCBI Gene Entry: [613](#)

UniProt Entry: [P11274](#)

Application Information

Molecular Weight: Predicted, 142 kDa; observed, 142,160 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat, monkey

Applications Tested: Western blotting (WB), immunohistochemistry (IHC), immunocytochemistry (IC)

Immunogen

A synthesized peptide derived from human Phospho-BCR (Y360)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Immunohistochemistry (IHC): 1:100-1:200

Immunocytochemistry (IC): 1:100-1:500

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

