

Anti-Phospho-CD115 (Y809) Rabbit Polyclonal Antibody



Catalog #: U0620

Aliases

FMS; Macrophage colony-stimulating factor 1 receptor; CSF-1 receptor; CSF-1-R; CSF-1R; M-CSF-R; Proto-oncogene c-Fms; CD115

Background

Gene Name: CSF1R

NCBI Gene Entry: [1436](#)

UniProt Entry: [P07333](#)

Application Information

Molecular Weight: Predicted, 107 kDa; observed, 140 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, bovine, dog, pig, sheep

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Phospho-CD115 (Y809)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of CD115 (Phospho-Y809) expression in MCF7 (A), HEK293T (B) whole cell lysates.