

Anti-Phospho-Beta-catenin (T41/S45) Rabbit Polyclonal Antibody



Catalog #: U0622

Aliases

CTNNB; Catenin beta-1; Beta-catenin

Background

Gene Name: CTNNB1

NCBI Gene Entry: [1499](#)

UniProt Entry: [P35222](#)

Application Information

Molecular Weight: Predicted, 85 kDa; observed, 92 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat, bovine, dog, pig, rabbit, sheep, zebrafish

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Phospho-Beta-catenin (T41/S45)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-Phospho-Beta-catenin (T41/S45) Rabbit Polyclonal Antibody

PAGE 2



Western blotting analysis of Beta-catenin (Phospho-T41/S45) expression in HeLa (A), HepG2 (B) whole cell lysates.

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM