

Anti-Phospho-LAT (Y200) Rabbit Polyclonal Antibody



Catalog #: U0651

Aliases

Linker for activation of T-cells family member 1; 36 kDa phospho-tyrosine adapter protein; pp36; p36-38

Background

Gene Name: LAT

NCBI Gene Entry: [27040](#)

UniProt Entry: [O43561](#)

Application Information

Molecular Weight: Predicted, 27 kDa; observed, 36, 38 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Phospho-LAT (Y200)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of LAT (Phospho-Y200) expression in Jurkat UV-treated (A), HeLa UV-treated (B) whole cell lysates.