

# Anti-Phospho-Progesterone Receptor (S190) Rabbit Polyclonal Antibody



**Catalog #: U0673**

## Aliases

NR3C3; Progesterone receptor; PR; Nuclear receptor subfamily 3 group C member 3

## Background

Gene Name: PGR

NCBI Gene Entry: [5241](#)

UniProt Entry: [P06401](#)

## Application Information

Molecular Weight: Predicted, 98 kDa; observed, 90, 118 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat, monkey

Applications Tested: Western blotting (WB), immunohistochemistry (IHC)

## Immunogen

A synthesized peptide derived from human Phospho-Progesterone Receptor (S190)

## Isotype

Rabbit IgG

## Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

## Storage

Store at -20 °C for one year.

## Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Immunohistochemistry (IHC): 1:50-1:100

**Note:** This product is for research use only.

## Validation Data

---

### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041


### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041


[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)

## Anti-Phospho-Progesterone Receptor (S190) Rabbit Polyclonal Antibody

PAGE 2



Immunohistochemical analysis of Progesterone Receptor (Phospho-S190) staining in human breast cancer formalin fixed paraffin embedded tissue section. The section was pre-treated using heat mediated antigen retrieval with sodium citrate buffer (pH 6.0). The section was then incubated with the antibody at room temperature and detected using an HRP conjugated compact polymer system. DAB was used as the chromogen. The section was then counterstained with haematoxylin and mounted with DPX.



Western blotting analysis of Progesterone Receptor (Phospho-S190) expression in mouse muscle (A), mouse lung (B), rat muscle (C) whole cell lysates.

---

### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)