

Anti-Phospho-PKC alpha (Y658) Rabbit Polyclonal Antibody



Catalog #: U0677

Aliases

PKCA; PRKACA; Protein kinase C alpha type; PKC-A; PKC-alpha

Background

Gene Name: PRKCA

NCBI Gene Entry: [5578](#)

UniProt Entry: [P17252](#)

Application Information

Molecular Weight: Predicted, 76 kDa; observed, 80 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat, bovine

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Phospho-PKC alpha (Y658)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-Phospho-PKC alpha (Y658) Rabbit Polyclonal Antibody

PAGE 2



Western blotting analysis of PKC alpha (Phospho-Y658) expression in HEK293T (A), HepG2 (B) whole cell lysates.

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM