

# Anti-Phospho-TIE2 (Y992) Rabbit Polyclonal Antibody



**Catalog #: U0695**

## Aliases

TIE2; VMCM; VMCM1; Angiopoietin-1 receptor; Endothelial tyrosine kinase; Tunica interna endothelial cell kinase; Tyrosine kinase with Ig and EGF homology domains-2; Tyrosine-protein kinase receptor TEK; Tyrosine-protein kinase receptor TIE-2; hTIE2; p140 TEK; CD202b

## Background

Gene Name: TEK

NCBI Gene Entry: [7010](#)

UniProt Entry: [Q02763](#)

## Application Information

Molecular Weight: Predicted, 125 kDa; observed, 160 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, bovine, zebrafish

Applications Tested: Western blotting (WB)

## Immunogen

A synthesized peptide derived from human Phospho-TIE2 (Y992)

## Isotype

Rabbit IgG

## Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

## Storage

Store at -20 °C for one year.

## Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

**Note:** This product is for research use only.

## Validation Data

---


### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)



Western blotting analysis of TIE2 (Phospho-Y992) expression in HuvEc PMA-treated (A) whole cell lysates.