

Anti-Phospho-ATF1 (S63) Rabbit Polyclonal Antibody



Catalog #: U0721

Aliases

Cyclic AMP-dependent transcription factor ATF-1; cAMP-dependent transcription factor ATF-1; Activating transcription factor 1; Protein TREB36

Background

Gene Name: ATF1

NCBI Gene Entry: [466](#)

UniProt Entry: [P18846](#)

Application Information

Molecular Weight: Predicted, 29 kDa; observed, 35 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, bovine

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Phospho-ATF1 (S63)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of ATF1 (Phospho-S63) expression in HEK293T (A), NIH3T3 (B), mouse kidney (C) whole cell lysates.