

Anti-Phospho-Tensin 2 (Y483) Rabbit Polyclonal Antibody



Catalog #: U0728

Aliases

KIAA1075; TNS2; Tensin-like C1 domain-containing phosphatase; C1 domain-containing phosphatase and tensin homolog; C1-TEN; Tensin-2

Background

Gene Name: TENC1

NCBI Gene Entry: [23371](#)

UniProt Entry: [Q63HR2](#)

Application Information

Molecular Weight: Predicted, 152 kDa; observed, 152 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Phospho-Tensin 2 (Y483)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-Phospho-Tensin 2 (Y483) Rabbit Polyclonal Antibody

PAGE 2



Western blotting analysis of Tensin 2 (Phospho-Y483) expression in A549 (A), HeLa (B), NIH3T3 (C) whole cell lysates.

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM