

Anti-Phospho-COP1 (S387) Rabbit Polyclonal Antibody



Catalog #: U0730

Aliases

COP1; RNF200; E3 ubiquitin-protein ligase RFWD2; Constitutive photomorphogenesis protein 1 homolog; hCOP1; RING finger and WD repeat domain protein 2; RING finger protein 200

Background

Gene Name: RFWD2

NCBI Gene Entry: [64326](#)

UniProt Entry: [Q8NHY2](#)

Application Information

Molecular Weight: Predicted, 80 kDa; observed, 110 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Phospho-COP1 (S387)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of COP1 (Phospho-S387) expression in Jurkat UV-treated (A), 293 UV-treated (B), K562 UV-treated (C) whole cell lysates.