

Anti-Phospho-MERTK/TYRO3 (Y753/685) Rabbit Polyclonal Antibody



Catalog #: U0753

Aliases

MERTK; MER; Tyrosine-protein kinase Mer; Proto-oncogene c-Mer; Receptor tyrosine kinase MerTK; TYRO3; BYK; DTK; RSE; SKY; Tyrosine-protein kinase receptor TYRO3; Tyrosine-protein kinase DTK; Tyrosine-protein kinase RSE; Tyrosine-protein kinase SKY; Tyrosine-protein kinase byk

Background

Gene Name: MERTK/TYRO3

NCBI Gene Entry: [10461/7301](#)

UniProt Entry: [Q12866/Q06418](#)

Application Information

Molecular Weight: Predicted, 110, 96 kDa; observed, 96 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat, chicken, zebrafish

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Phospho-MERTK/TYRO3 (Y753/685)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:2,000

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-Phospho-MERTK/TYRO3 (Y753/685) Rabbit Polyclonal Antibody

PAGE 2



Western blotting analysis of MERTK/TYRO3 (Phospho-Y753/685) expression in HEK293T (A), H1688 (B) whole cell lysates.

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM