

Anti-Phospho-BLK (Y389) Rabbit Polyclonal Antibody



Catalog #: U0756

Aliases

Tyrosine-protein kinase Blk; B lymphocyte kinase; p55-Blk

Background

Gene Name: BLK

NCBI Gene Entry: [640](#)

UniProt Entry: [P51451](#)

Application Information

Molecular Weight: Predicted, 57 kDa; observed, 58 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Phospho-BLK (Y389)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of BLK (Phospho-Y389) expression in HEK293T (A), PANC1 (B), CT26 (C), rat spleen (D) whole cell lysates.