

Anti-Phospho-CD45 (Y1216) Rabbit Polyclonal Antibody



Catalog #: U0761

Aliases

CD45; Receptor-type tyrosine-protein phosphatase C; Leukocyte common antigen; L-CA; T200; CD45

Background

Gene Name: PTPRC

NCBI Gene Entry: [5788](#)

UniProt Entry: [P08575](#)

Application Information

Molecular Weight: Predicted, 147 kDa; observed, 180 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Phospho-CD45 (Y1216)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of CD45 (Phospho-Y1216) expression in HEK293T (A), MCF7 (B), PC12 (C), AML12 (D) whole cell lysates.