

Anti-Phospho-PKC beta (T641) Rabbit Polyclonal Antibody



Catalog #: U0779

Aliases

PRKCB; PKCB; PRKCB1; Protein kinase C beta type; PKC-B; PKC-beta

Background

Gene Name: PRKCB

NCBI Gene Entry: [5579](#)

UniProt Entry: [P05771](#)

Application Information

Molecular Weight: Predicted, 76 kDa; observed, 80 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat, bovine, zebrafish

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Phospho-PKC beta (T641)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-Phospho-PKC beta (T641) Rabbit Polyclonal Antibody

PAGE 2



Western blotting analysis of PKC beta (Phospho-T641) expression in Hela (A), HEK293T (B), mouse brain (C), rat brain (D) whole cell lysates.

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM