

Anti-Phospho-eNOS (S615) Rabbit Polyclonal Antibody



Catalog #: U0804

Aliases

Nitric oxide synthase endothelial; Constitutive NOS; cNOS; EC-NOS; Endothelial NOS; eNOS; NOS type III; NOSIII

Background

Gene Name: NOS3

NCBI Gene Entry: [4846](#)

UniProt Entry: [P29474](#)

Application Information

Molecular Weight: Predicted, 133, 69 kDa; observed, 140 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat, bovine, dog, pig

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Phospho-eNOS (S615)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of eNOS (Phospho-S615) expression in K562 (A), mouse spleen (B), rat spleen (C) whole cell lysates.