

Anti-Phospho-Ubiquitin (S65) Rabbit Polyclonal Antibody



Catalog #: U0807

Aliases

UBB; Polyubiquitin-B; UBC; Polyubiquitin-C; RPS27A; UBA80; UBCEP1; Ubiquitin-40S ribosomal protein S27a; Ubiquitin carboxyl extension protein 80; UBA52; UBCEP2; Ubiquitin-60S ribosomal protein L40; CEP52; Ubiquitin A-52 residue ribosomal protein fusion product 1

Background

Gene Name: UBA52/RPS27A/UBB/UBC

NCBI Gene Entry: [7314/7316/6233/7311](#)

UniProt Entry: [P0CG47/P0CG48/P62979/P62987](#)

Application Information

Molecular Weight: Predicted, 25, 77, 17, 14 kDa; observed, 15 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Phospho-Ubiquitin (S65)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of Ubiquitin (Phospho-S65) expression in C6 (A), AML12 (B), SGC7901-UV treated (C), SGC7901 (D) whole cell lysates.