

Anti-Phospho-PINK1 (S228) Rabbit Polyclonal Antibody



Catalog #: U0808

Aliases

Serine/threonine-protein kinase PINK1 mitochondrial; BRPK; PTEN-induced putative kinase protein 1

Background

Gene Name: PINK1

NCBI Gene Entry: [65018](#)

UniProt Entry: [Q9BXM7](#)

Application Information

Molecular Weight: Predicted, 62 kDa; observed, 60 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB), immunohistochemistry (IHC)

Immunogen

A synthesized peptide derived from human Phospho-PINK1 (S228)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Immunohistochemistry (IHC): 1:50-1:100

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of PINK1 (Phospho-S228) expression in mouse testis (A), rat testis (B), HEK293T (C), MCF7 (D), U87MG (E) whole cell lysates.