

Anti-Phospho-Histone Deacetylase 2 (S394) Rabbit Polyclonal Antibody



Catalog #: U0810

Aliases

Histone deacetylase 2; HD2

Background

Gene Name: HDAC2

NCBI Gene Entry: [3066](#)

UniProt Entry: [Q92769](#)

Application Information

Molecular Weight: Predicted, 55 kDa; observed, 60 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat, chicken

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Phospho-Histone Deacetylase 2 (S394)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-Phospho-Histone Deacetylase 2 (S394) Rabbit Polyclonal Antibody

PAGE 2



Western blotting analysis of Histone Deacetylase 2 (Phospho-S394) expression in HEK293T (A), H1792 (B), Panc1 (C) whole cell lysates.

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM