

Anti-Phospho-FRS2 (Y196) Rabbit Polyclonal Antibody



Catalog #: U0815

Aliases

Fibroblast growth factor receptor substrate 2; FGFR substrate 2; FGFR-signaling adaptor SNT; Suc1-associated neurotrophic factor target 1; SNT-1

Background

Gene Name: FRS2

NCBI Gene Entry: [10818](#)

UniProt Entry: [Q8WU20](#)

Application Information

Molecular Weight: Predicted, 57 kDa; observed, 57, 75 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Phospho-FRS2 (Y196)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of FRS2 (Phospho-Y196) expression in HepG2 (A), HEK293T (B) whole cell lysates.