

# Anti-Phospho-cTnI (S23/S24) Rabbit Polyclonal Antibody



**Catalog #: U0828**

## Aliases

TNNC1; Troponin I cardiac muscle; Cardiac troponin I

## Background

Gene Name: TNNI3

NCBI Gene Entry: [7137](#)

UniProt Entry: [P19429](#)

## Application Information

Molecular Weight: Predicted, 24 kDa; observed, 24 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

## Immunogen

A synthesized peptide derived from human Phospho-cTnI (S23/S24)

## Isotype

Rabbit IgG

## Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

## Storage

Store at -20 °C for one year.

## Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

**Note:** This product is for research use only.

## Validation Data

---


### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)



Western blotting analysis of cTnI (Phospho-S23/S24) expression in rat heart (A), H9C2 (B), A549 (C), HEK293T (D), mouse testis (E) whole cell lysates.