

Anti-Phospho-E Cadherin (S844) Rabbit Polyclonal Antibody



Catalog #: U0829

Aliases

CDHE; UVO; Cadherin-1; CAM 120/80; Epithelial cadherin; E-cadherin; Uvomorulin; CD324

Background

Gene Name: CDH1

NCBI Gene Entry: [999](#)

UniProt Entry: [P12830](#)

Application Information

Molecular Weight: Predicted, 97 kDa; observed, 120 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat, dog

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Phospho-E Cadherin (S844)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of E Cadherin (Phospho-S844) expression in HEK293T (A), HepG2 (B), A549 (C), PC12 (D) whole cell lysates.