

Anti-Phospho-MKK5 (S311) Rabbit Polyclonal Antibody



Catalog #: U0830

Aliases

MEK5; MKK5; PRKMK5; Dual specificity mitogen-activated protein kinase kinase 5; MAP kinase kinase 5; MAPKK 5; MAPK/ERK kinase 5; MEK 5

Background

Gene Name: MAP2K5

NCBI Gene Entry: [5607](#)

UniProt Entry: [Q13163](#)

Application Information

Molecular Weight: Predicted, 50 kDa; observed, 50, 46 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Phospho-MKK5 (S311)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of MKK5 (Phospho-S311) expression in mouse heart (A), rat heart (B), H1792 (C), MCF7 (D) whole cell lysates.