

Anti-Phospho-NF-kappaB p105 (S903) Rabbit Polyclonal Antibody



Catalog #: U0837

Aliases

Nuclear factor NF-kappa-B p105 subunit; DNA-binding factor KBF1; EBP-1; Nuclear factor of kappa light polypeptide gene enhancer in B-cells 1

Background

Gene Name: NFKB1

NCBI Gene Entry: [4790](#)

UniProt Entry: [P19838](#)

Application Information

Molecular Weight: Predicted, 105 kDa; observed, 120 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, rat

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Phospho-NF-kappaB p105 (S903)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of NF-kappaB p105 (Phospho-S903) expression in MCF7-TNFa (A), rat spleen (B) whole cell lysates.