

Anti-Phospho-Progesterone Receptor (S400) Rabbit Polyclonal Antibody



Catalog #: U0838

Aliases

NR3C3; Progesterone receptor; PR; Nuclear receptor subfamily 3 group C member 3

Background

Gene Name: PGR

NCBI Gene Entry: [5241](#)

UniProt Entry: [P06401](#)

Application Information

Molecular Weight: Predicted, 98 kDa; observed, 118 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB), immunocytochemistry (IC)

Immunogen

A synthesized peptide derived from human Phospho-Progesterone Receptor (S400)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Immunocytochemistry (IC): 1:100-1:500

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041


WWW.GENUINBIOTECH.COM

Anti-Phospho-Progesterone Receptor (S400) Rabbit Polyclonal Antibody

PAGE 2



Immunocytochemical analysis of Progesterone Receptor (Phospho-S400) staining in A549 cells. Formalin-fixed cells were permeabilized with 0.1% Triton X-100 in TBS for 5-10 minutes and blocked with 3% BSA-PBS for 30 minutes at room temperature. Cells were probed with the primary antibody in 3% BSA-PBS and incubated overnight at 4 °C in a humidified chamber. Cells were washed with PBST and incubated with a Alexa Fluor 594-conjugated secondary antibody (red) in PBS at room temperature in the dark.



Western blotting analysis of Progesterone Receptor (Phospho-S400) expression in mouse testis (A), rat ovary (B), A2780 (C), HeLa (D) whole cell lysates.

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM