

Anti-Phospho-FOXO3 (S315) Rabbit Polyclonal Antibody



Catalog #: U0875

Aliases

FKHRL1; FOXO3A; Forkhead box protein O3; AF6q21 protein; Forkhead in rhabdomyosarcoma-like 1

Background

Gene Name: FOXO3

NCBI Gene Entry: [2309](#)

UniProt Entry: [O43524](#)

Application Information

Molecular Weight: Predicted, 71 kDa; observed, 82, 97 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Phospho-FOXO3 (S315)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of FOXO3 (Phospho-S315) expression in HeLa (A), NIH3T3 (B) whole cell lysates.