

Anti-Phospho-USP28 (S714) Rabbit Polyclonal Antibody



Catalog #: U0876

Aliases

KIAA1515; Ubiquitin carboxyl-terminal hydrolase 28; Deubiquitinating enzyme 28; Ubiquitin thioesterase 28; Ubiquitin-specific-processing protease 28

Background

Gene Name: USP28

NCBI Gene Entry: [57646](#)

UniProt Entry: [Q96RU2](#)

Application Information

Molecular Weight: Predicted, 122 kDa; observed, 122 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Phospho-USP28 (S714)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of USP28 (Phospho-S714) expression in HeLa (A), A549 (B) whole cell lysates.