

Anti-Phospho-ZAP70 (Y292) Rabbit Polyclonal Antibody



Catalog #: U0884

Aliases

SRK; Tyrosine-protein kinase ZAP-70; 70 kDa zeta-chain associated protein; Syk-related tyrosine kinase

Background

Gene Name: ZAP70

NCBI Gene Entry: [7535](#)

UniProt Entry: [P43403](#)

Application Information

Molecular Weight: Predicted, 69 kDa; observed, 70 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Phospho-ZAP70 (Y292)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of ZAP70 (Phospho-Y292) expression in Jurkat (A), Myla2059 (B), A549 (C), rat thymus (D) whole cell lysates.