

Anti-Phospho-NCOA1 (T1179) Rabbit Polyclonal Antibody



Catalog #: U0909

Aliases

BHLHE74; SRC1; Nuclear receptor coactivator 1; NCoA-1; Class E basic helix-loop-helix protein 74; bHLHe74; Protein Hin-2; RIP160; Renal carcinoma antigen NY-REN-52; Steroid receptor coactivator 1; SRC-1

Background

Gene Name: NCOA1

NCBI Gene Entry: [8648](#)

UniProt Entry: [Q15788](#)

Application Information

Molecular Weight: Predicted, 156 kDa; observed, 180 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat, pig

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Phospho-NCOA1 (T1179)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-Phospho-NCOA1 (T1179) Rabbit Polyclonal Antibody

PAGE 2



Western blotting analysis of NCOA1 (Phospho-T1179) expression in CT26 (A), U87MG (B) whole cell lysates.

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM