

Anti-Phospho-TRK B (Y817) Rabbit Polyclonal Antibody



Catalog #: U0918

Aliases

TRKB; BDNF/NT-3 growth factors receptor; GP145-TrkB; Trk-B; Neurotrophic tyrosine kinase receptor type 2; TrkB tyrosine kinase; Tropomyosin-related kinase B

Background

Gene Name: NTRK2

NCBI Gene Entry: [4915](#)

UniProt Entry: [Q16620](#)

Application Information

Molecular Weight: Predicted, 91 kDa; observed, 92 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat, chicken

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Phospho-TRK B (Y817)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of TRK B (Phospho-Y817) expression in U87MG (A), HEK293T (B) whole cell lysates.