

Anti-Phospho-DAB1 (Y220) Rabbit Polyclonal Antibody



Catalog #: U0926

Aliases

Disabled homolog 1

Background

Gene Name: DAB1

NCBI Gene Entry: [1600](#)

UniProt Entry: [O75553](#)

Application Information

Molecular Weight: Predicted, 63 kDa; observed, 70 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Phospho-DAB1 (Y220)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of DAB1 (Phospho-Y220) expression in mouse brain (A), rat brain (B) whole cell lysates.