

Anti-Phospho-Topoisomerase 2 alpha (S1213) Rabbit Polyclonal Antibody



Catalog #: U0933

Aliases

TOP2; DNA topoisomerase 2-alpha; DNA topoisomerase II alpha isozyme

Background

Gene Name: TOP2A

NCBI Gene Entry: [7153](#)

UniProt Entry: [P11388](#)

Application Information

Molecular Weight: Predicted, 174 kDa; observed, 190 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Phospho-Topoisomerase 2 alpha (S1213)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-Phospho-Topoisomerase 2 alpha (S1213) Rabbit Polyclonal Antibody

PAGE 2



Western blotting analysis of Topoisomerase 2 alpha (Phospho-S1213) expression in SHSY5Y (A), Myla2059 (B) whole cell lysates.

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM