

# Anti-Phospho-14-3-3 sigma (S186) Rabbit Polyclonal Antibody



**Catalog #: U0937**

## Aliases

HME1; 14-3-3 protein sigma; Epithelial cell marker protein 1; Stratifin

## Background

Gene Name: SFN

NCBI Gene Entry: [2810](#)

UniProt Entry: [P31947](#)

## Application Information

Molecular Weight: Predicted, 27 kDa; observed, 30 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat, bovine, pig, sheep

Applications Tested: Western blotting (WB)

## Immunogen

A synthesized peptide derived from human Phospho-14-3-3 sigma (S186)

## Isotype

Rabbit IgG

## Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

## Storage

Store at -20 °C for one year.

## Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

**Note:** This product is for research use only.

## Validation Data

---

### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041


### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)

## Anti-Phospho-14-3-3 sigma (S186) Rabbit Polyclonal Antibody

PAGE 2



Western blotting analysis of 14-3-3 sigma (Phospho-S186) expression in COS7 (A), MEF (B), SHSY5Y (C), HCT116 (D), A549 (E) whole cell lysates.

### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)