

Anti-Phospho-PKC gamma (T655) Rabbit Polyclonal Antibody



Catalog #: U0956

Aliases

PKCG; Protein kinase C gamma type; PKC-gamma

Background

Gene Name: PRKCG

NCBI Gene Entry: [5582](#)

UniProt Entry: [P05129](#)

Application Information

Molecular Weight: Predicted, 78 kDa; observed, 80 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Phospho-PKC gamma (T655)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of PKC gamma (Phospho-T655) expression in mouse brain (A), rat brain (B) whole cell lysates.