

Anti-Phospho-Presenilin 1 (S353) Rabbit Polyclonal Antibody



Catalog #: U0966

Aliases

AD3; PS1; PSNL1; Presenilin-1; PS-1; Protein S182

Background

Gene Name: PSEN1

NCBI Gene Entry: [5663](#)

UniProt Entry: [P49768](#)

Application Information

Molecular Weight: Predicted, 52 kDa; observed, 52, 22 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Phospho-Presenilin 1 (S353)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-Phospho-Presenilin 1 (S353) Rabbit Polyclonal Antibody

PAGE 2



Western blotting analysis of Presenilin 1 (Phospho-S353) expression in A549 (A), THP1 (B) whole cell lysates.

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM