

Anti-Phospho-Von Hippel Lindau (S68) Rabbit Polyclonal Antibody



Catalog #: U0967

Aliases

Von Hippel-Lindau disease tumor suppressor; Protein G7; pVHL

Background

Gene Name: VHL

NCBI Gene Entry: [7428](#)

UniProt Entry: [P40337](#)

Application Information

Molecular Weight: Predicted, 24 kDa; observed, 24 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Phospho-Von Hippel Lindau (S68)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of Von Hippel Lindau (Phospho-S68) expression in HEK293T (A), MCF7 (B), HeLa (C) whole cell lysates.