

Anti-Phospho-Zyxin (S142) Rabbit Polyclonal Antibody



Catalog #: U0968

Aliases

Zyxin; Zyxin-2

Background

Gene Name: ZYX

NCBI Gene Entry: [7791](#)

UniProt Entry: [Q15942](#)

Application Information

Molecular Weight: Predicted, 61 kDa; observed, 61 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Phospho-Zyxin (S142)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of Zyxin (Phospho-S142) expression in HuT78 (A), mouse muscle (B), rat lung (C), rat muscle (D) whole cell lysates.