

Anti-Phospho-MYL2 (S15) Rabbit Polyclonal Antibody



Catalog #: U0970

Aliases

Myosin regulatory light chain 2 ventricular/cardiac muscle isoform; MLC-2; MLC-2v

Background

Gene Name: MYL2

NCBI Gene Entry: [4633](#)

UniProt Entry: [P10916](#)

Application Information

Molecular Weight: Predicted, 18 kDa; observed, 18 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, rat

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Phospho-MYL2 (S15)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of MYL2 (Phospho-S15) expression in HEK293T (A), PC12 (B), SGC7901 (C) whole cell lysates.