

Anti-Phospho-STAT3 (S727) Rabbit Monoclonal Antibody



Catalog #: U1056

Aliases

APRF; Signal transducer and activator of transcription 3; Acute-phase response factor

Background

Gene Name: STAT3

NCBI Gene Entry: [6774](#)

UniProt Entry: [P40763](#)

Application Information

Molecular Weight: Predicted, 88 kDa; observed, 88 kDa

Clonality: Rabbit monoclonal antibody

Clone ID: 24GB7820

Species Reactivity: Human, rat

Applications Tested: Western blotting (WB), immunocytochemistry (IC)

Immunogen

A synthesized peptide derived from human Phospho-STAT3 (S727)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, 0.05% BSA and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Immunocytochemistry (IC): 1:50-1:100

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of STAT3 (Phospho-S727) expression in Jurkat (A), C6 (B), HeLa (C) whole cell lysates.