

Anti-Phospho-JNK1/2/3 (T183/Y185) Rabbit Monoclonal Antibody



Catalog #: U1060

Aliases

JNK1; PRKM8; SAPK1; SAPK1C; Mitogen-activated protein kinase 8; MAP kinase 8; MAPK 8; JNK-46; Stress-activated protein kinase 1c; SAPK1c; Stress-activated protein kinase JNK1; c-Jun N-terminal kinase 1

Background

Gene Name: MAPK8

NCBI Gene Entry: [5599](#)

UniProt Entry: [P45983](#)

Application Information

Molecular Weight: Predicted, 48 kDa; observed, 46,54 kDa

Clonality: Rabbit monoclonal antibody

Clone ID: 24GB7840

Species Reactivity: Human, rat

Applications Tested: Western blotting (WB), immunoprecipitation (IP)

Immunogen

A synthesized peptide derived from human Phospho-JNK1/2/3 (T183/Y185)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, 0.05% BSA and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Immunoprecipitation (IP): 1:10-1:50

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of JNK1/2/3 (Phospho-T183/Y185) expression in HeLa (A), A549 (B), C6 (C) whole cell lysates.