

# Anti-Phospho-c-Met (Y1234/Y1235) Mouse Monoclonal Antibody



**Catalog #: U1065**

## Aliases

Hepatocyte growth factor receptor; HGF receptor; HGF/SF receptor; Proto-oncogene c-Met; Scatter factor receptor; SF receptor; Tyrosine-protein kinase Met

## Background

Gene Name: MET

NCBI Gene Entry: [4233](#)

UniProt Entry: [P08581](#)

## Application Information

Molecular Weight: Predicted, 155 kDa; observed, 145 kDa

Clonality: Mouse monoclonal antibody

Clone ID: 24GB7730

Species Reactivity: Human

Applications Tested: Western blotting (WB)

## Immunogen

A synthesized peptide derived from human Phospho-c-Met (Y1234/Y1235)

## Isotype

Mouse IgG1

## Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

## Storage

Store at -20 °C for one year.

## Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

**Note:** This product is for research use only.

## Validation Data

---

### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041


### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)

## Anti-Phospho-c-Met (Y1234/Y1235) Mouse Monoclonal Antibody

PAGE 2



Western blotting analysis of c-Met (Phospho-Y1234/Y1235) expression in HepG2 (A) whole cell lysates.

---

### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)