

Anti-Phospho-PLM (S68) Mouse Monoclonal Antibody



Catalog #: U1068

Aliases

PLM; Phospholemman; FXYP domain-containing ion transport regulator 1

Background

Gene Name: FXYP1

NCBI Gene Entry: [5348](#)

UniProt Entry: [O00168](#)

Application Information

Molecular Weight: Predicted, 10 kDa; observed, 8 kDa

Clonality: Mouse monoclonal antibody

Clone ID: 24GB7745

Species Reactivity: Mouse

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Phospho-PLM (S68)

Isotype

Mouse IgG1

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of PLM (Phospho-S68) expression in 293 transfected with a plasmid encoding PLM (A) whole cell lysates.