

Anti-Phospho-Erbin (Y1104) Rabbit Polyclonal Antibody



Catalog #: U1081

Aliases

ERBB2IP; ERBIN; KIAA1225; LAP2; Protein LAP2; Densin-180-like protein; ErbB2-interacting protein; Erbin

Background

Gene Name: ERBB2IP

NCBI Gene Entry: [55914](#)

UniProt Entry: [Q96RT1](#)

Application Information

Molecular Weight: Predicted, 158 kDa; observed, 172 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Phospho-Erbin (Tyr1104)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:2,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of lysates from HeLa cells, using phospho-Erbin (Phospho-Tyr1104) antibody.