

Anti-Phospho-Ob-R (Y1141) Rabbit Polyclonal Antibody



Catalog #: U1757

Aliases

LEPR; DB; OBR; Leptin receptor; LEP-R; HuB219; OB receptor; OB-R; CD antigen CD295

Background

Gene Name: LEPR/OBR

NCBI Gene Entry: [3953](#)

UniProt Entry: [P48357](#)

Application Information

Molecular Weight: Predicted, 132 kDa; observed, 132 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Phospho-Ob-R (Tyr1141)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:2,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of extracts from K562 cells, using Phospho-Ob-R (Y1141) Polyclonal Antibody.