

Anti-Phospho-ATP1A1 (Y260) Rabbit Polyclonal Antibody



Catalog #: U1768

Aliases

Sodium/potassium-transporting ATPase subunit alpha-1 (Na(+)/K(+) ATPase alpha-1 subunit) (EC 3.6.3.9) (Sodium pump subunit alpha-1)

Background

Gene Name: ATP1A1

NCBI Gene Entry: [476](#)

UniProt Entry: [P05023](#)

Application Information

Molecular Weight: Predicted, 113 kDa; observed, 115 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB), immunohistochemistry (IHC)

Immunogen

A synthesized peptide derived from human Phospho-Sodium Potassium ATPase alpha-1 (Tyr260)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:2,000

Immunohistochemistry-Paraffin (IHC-P): 1:50-1:300

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of lysates from 293 cells treated with PMA, using phospho-Na⁺/K⁺-ATPase α1 (Phospho-Tyr260) antibody.