

Anti-Phospho-Catenin delta-1 (Y904) Rabbit Polyclonal Antibody



Catalog #: U1820

Aliases

Catenin delta-1 (Cadherin-associated Src substrate) (CAS) (p120 catenin) (p120(ctn)) (p120(cas))

Background

Gene Name: CTNND1/KIAA0384

NCBI Gene Entry: [1500](#)

UniProt Entry: [O60716](#)

Application Information

Molecular Weight: Predicted, 108 kDa; observed, 108 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Phospho-Catenin δ -1 (Tyr904)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:1,000-1:2,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of HEK293, mouse brain, using primary antibody at 1:1,000 dilution.
Secondary antibody was diluted at 1:10,000