

Anti-Phospho-CYLD (S418) Rabbit Polyclonal Antibody



Catalog #: U1835

Aliases

Ubiquitin carboxyl-terminal hydrolase CYLD (EC 3.4.19.12) (Deubiquitinating enzyme CYLD) (Ubiquitin thioesterase CYLD) (Ubiquitin-specific-processing protease CYLD)

Background

Gene Name: CYLD/CYLD1/KIAA0849/HSPC057

NCBI Gene Entry: [1540](#)

UniProt Entry: [Q9NQC7](#)

Application Information

Molecular Weight: Predicted, 107 kDa; observed, 105 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Phospho-CYLD (Ser418)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:1,000-1:2,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of 1 A431 treated with LPS, 2 A431, using primary antibody at 1:1,000 dilution. Secondary antibody was diluted at 1:10,000