

Anti-Phospho-FLT3 (Y589/591) Rabbit Polyclonal Antibody



Catalog #: U1857

Aliases

Receptor-type tyrosine-protein kinase FLT3 (EC 2.7.10.1) (FL cytokine receptor) (Fetal liver kinase-2) (FLK-2) (Fms-like tyrosine kinase 3) (FLT-3) (Stem cell tyrosine kinase 1) (STK-1) (CD antigen CD135)

Background

Gene Name: FLT3/CD135/FLK2/STK1

NCBI Gene Entry: [2322](#)

UniProt Entry: [P36888](#)

Application Information

Molecular Weight: Predicted, 113 kDa; observed, 117 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Phospho-FLT3 (Tyr589/591)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:1,000-1:2,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of mouse brain , rat brain ,using primary antibody at 1:1,000 dilution.
Secondary antibody was diluted at 1:10,000